

## Methadone Cross-Reactivity Table

*The following compounds tested POSITIVE using the Methadone assay at the 300 ng/mL cutoff.*

Postive Compounds	Trade Name	Concentration Tested (ng/mL)
$\alpha$ -Methadol		33,333
LAAM		20,000
Methadol		25,000
Methadone	Dolophine	300

The following compounds tested NEGATIVE on the Methadone 300 ng/mL assay.

Negative Compounds	Trade Name	Concentration Tested (ng/mL)
1,3-Dimethylbarbituric acid		500,000
10,11-Dihydrocarbamazepine		500,000
1-Phenylcyclohexylamine		500,000
4-OH-PCP HCl		500,000
5,5-Diphenylhydantoin (Phenytoin)	Dilantin	500,000
5-Hydroxyindole-2-Carboxylic acid		500,000
5-Hydroxyindole-3-Acetic acid		500,000
5-Hydroxytryptamine (Serotonin)		500,000
6-Monoacetyl morphine		500,000
8,11-di-OH- $\Delta^9$ THC		100,000
11-nor- $\Delta^9$ -THC-COOH		10,000
Acetaminophen (paracetamol)	Anacin, Datri Extra, Liquiprin, Panadol, Tempra, Tylenol	500,000
Acetanilide		500,000
Acetylcysteine	Mucomyst, Mucosil, Parvalex	500,000
Acetazolamide	Diamox, Acetazolam	500,000
Acetophenetidine		500,000
N-Acetylprocainamide		500,000
Acetylsalicylic acid	Aspirin, Bufferin	500,000
Acyclovir	Zovirax	500,000
Albuterol	Salbutamol, Proventil, Ventolin	500,000
Allobarbitol	Asmac, Cibalgine, Diadol, Dially Barbituric Acid, Dial, Malilum, Sornnocodal	500,000
Allopurinol	Alloprim, Apo-Allopurinol, Lopurin, Novopurol, Purhol	500,000
Alpha-Hydroxyalprazolam		500,000
Alphaprodine HCl		500,000
Alphenal	Alphenate, Allofenyl, Efrogal, Fenallymal	500,000
Alprazolam	Tafil, Valeans, Xanax, Xandor	500,000
Alprenolol HCl		500,000
Amantadine HCl	Symadine, Symmentrel	500,000
p-Aminobenzoic acid (PABA)		500,000
Aminoglutethimide	Cytadren, Orimeten	500,000
Aminopyrine		500,000

The following compounds tested NEGATIVE on the Methadone 300 ng/mL assay.

Negative Compounds	Trade Name	Concentration Tested (ng/mL)
Amitriptyline HCl	Elavil, Endep, Enovil, Emitrip	500,000
Ammonium chloride		500,000
Amobarbital	Amytal	500,000
Amoxapine	Asendin	500,000
Amoxicillin	Amoxil, Polymor, Trimox, Wymox	100,000
Amphetamine		500,000
Ampicillin Na <sup>+</sup> Salt	Omnipen, Paracetamol, Polycillin, Principen	500,000
Ampicillin anhydrous		500,000
Aniline HCl		500,000
Aprobarbital	Alurate	500,000
l-Ascorbic acid	Ascorbicap, Cebid, Cecon, Cevalin, Flavorcee, Nice, Vita-C	500,000
Aspartame	Equal, Nutrasweet	500,000
Astemizole	Hismanal	500,000
Atenolol	Tenormin	500,000
Atropine	Atropen, Attropisol, Atrosulf, Saltropine	500,000
Barbital	Deba	50,000
Barbituric acid		500,000
Beclomethasone dipropionate	Beclovent	500,000
Bemegride		500,000
Benactyzine		500,000
Benzilic acid		500,000
Benzocaine	Anbesol, Americaine, Otic, Hurrricane, Orajel, Zilactin	500,000
Benzoic acid	Atrosept, Dolsed, Prosed, UAA, Uridon, Urised, Uritab, Uritin	500,000
Benzonatate	Tessalon	500,000
Benzophenone		500,000
Benzotropine mesylate	Cogentin	500,000
Benzoylecgonine		500,000
Benzylamine HCl		500,000
Bethanechol chloride	Duvoid, Urabeth	500,000
Bromazepam	Lectopam	500,000
Bromocriptine mesylate	Ergoset, Parlodel	500,000
Bupivacaine	Marcaine, Sensorcaine	500,000
Buprenorphine	Buprenex	500,000

The following compounds tested NEGATIVE on the Methadone 300 ng/mL assay.

Negative Compounds	Trade Name	Concentration Tested (ng/mL)
Buspirone HCl	Buspar	500,000
Butabarbital	Butisol, Buticaps, Busodium	500,000
Butalbital	Butalgen, Fiorgen, Fiorinal, Fiormor, Isobutyl, Isolin, Isollyl, Laninoif	500,000
Butethal	Etoval	500,000
Butorphanol	Stadol	500,000
Butyrophenone		500,000
Caffeine	Coffe-Break, Durivitan, No-Doz, Percoffedrinol N, Pro-Plus, Vivarin	500,000
Calcium hypochlorite		500,000
Captopril	Capoten	500,000
Carbamazepine	Atretol, Carbatrel, Tegretol	500,000
Carbromal		500,000
Carisoprodol	Rela, Sodol, Soma, Soprodol, Soridel	500,000
Cefadroxil monohydrate	Duricet, Ultracef	500,000
Cefuroxime	Zinacef	500,000
Cephalexin	Cefanex, Keflex	500,000
Cephaloridine		500,000
Cephradine	Velosef, Anspor	500,000
Chloramphenicol	Chloroptic, Chloromycetin, Ocu-Chlor	500,000
Chlordiazepoxide	Libritabs, Librium, Sereen	100,000
Chloroquine	Aralen	500,000
Chlorothiazide	Diuril	500,000
Chlorotrianisene	Tace	500,000
Chlorpromazine HCl	Thorazine	500,000
Chlorpropamide	Diabinese, Glucamide	500,000
Chlorprothixene	Taractan	500,000
Chlorzoxazone	Paraflex, Parafon Forte DSC	500,000
Cholesterol		500,000
Cimetidine	Tagamet	500,000
Cinoxacin	Cinobac	500,000
Ciprofloxacin HCl	Cipnol	500,000
Clindamycin	Cleocin	500,000
Clomipramine HCl	Anafranil	500,000
Clonazepam	Klonopin	500,000
Clotrimazole	Mycelex	500,000

The following compounds tested NEGATIVE on the Methadone 300 ng/mL assay.

Negative Compounds	Trade Name	Concentration Tested (ng/mL)
Cloxacillin	Cloxapen, Dynapen, Pathocil, Tegopen	500,000
Cocaine		500,000
Codeine		500,000
Colchicine		500,000
Cortisol	A-HydroCort, Hydrocortone Acetate, Solu-Cortef, Hydrocortisone	500,000
Cortisone	Cortone	500,000
Cortolone		500,000
Cotinine		500,000
Creatine hemisulphate		500,000
Creatine phosphate		500,000
Cyclobenzaprine	Cycloflex, Flexeril	500,000
Cyclohexylamine		500,000
Cyclopentobarbital		500,000
Cyclophosphamide	Cytoxan, Neosar	500,000
Cyproheptadine	Periactin	500,000
Demoxepam		500,000
Deoxycorticosterone		500,000
Deoxyepinephrine HCl		500,000
Desalkylfurazepam		500,000
Desipramine HCl	Norpramin, Petrofrane	500,000
o-desmethylvenlafaxine		1,000,000
Dextromethorphan HBr	Benylin DM, Delsym, Hold, Perfussin 8 hr, Mediquell, Sucrets	500,000
Diacetyl morphine	Heroin	500,000
Diazepam	Valium, Valrelease, Vazepam	500,000
Dibucaine HCl	Nipercainal	500,000
Diclofenac	Cataflam, Pennsaid, Voltaren	500,000
Dicyclomine HCl	Antipas, Bentyl Hydrochloride, Dibent, Dilomine	500,000
Diethylsuccinate		500,000
Diflunisal	Dolobid	500,000
Digoxin	Lanoxin	100,000
Dihydrocarbamazepine		500,000
Dihydrocodeine	DHC Plus, SynalgospDC	100,000
Diltiazem HCl	Cardizem, Dilacor, Dilitazem, Tiazac	500,000

The following compounds tested NEGATIVE on the Methadone 300 ng/mL assay.

Negative Compounds	Trade Name	Concentration Tested (ng/mL)
Diphenhydramine HCl	Banophen, Benadryl, Diphedryl, Simply Sleep	500,000
Diphenoxylate HCl	Lomotil	500,000
Dipyridamole	Persantine, Permole	500,000
Disopyramide phosphate	Norpace, Rythmodan	12,000
Disulfiram	Antabuse	500,000
Dopa	Dopar, Larodopa	500,000
Dopamine HCl	Dopastat, Intropin	500,000
Doxepin HCl	Adapin, Sinequan	500,000
Doxylamine	Decapryn, Nighttime Sleep Aid, Unisom	500,000
Dronabinol ( $\Delta^9$ -THC)	Marinol	500,000
Ecgonine		500,000
Ecgonine methyl ester		500,000
EDDP (Primary Methadone Metabolite)		500,000
EMDP (Secondary Methadone Metabolite)		100,000
Enalapril	Vasotec	500,000
Ephedrine	Kondons Nasal, Pretz-D	500,000
Epinephrine	Adrenalin, AsthmaHaler, Bronitin, Epifrin, Primatene, Vaponefrin	500,000
Erythromycin	AK-Mycin, E-Mycin, Eryc, Ery-Tabs, Ilosone, Wyamycin	500,000
Esterone		500,000
$\beta$ -Estradiol (E-2)	Alora, Climara, Delestrogen, DepGynogen, Depo-Estradiol, Depogen, E-Cypionae, Esclim, Estinyl, Estrace, Estraderm, Estragyn LA 5, Estring, EstroCyp, Estro-L.A., FemPatch, Gynodiol, Gynodiol, Gynogen L.A., Innofem, Menaval, Noven, Vagifem, Valergen, Vivelle, Vievelle-Dot	500,000
Estriol (E-3)		500,000
Estrone-3-Sulfate, Na <sup>+</sup> (E-1)	Estragyn, Estratab, Kestrone, Menest, Ogen, Ortho-Est	500,000
Ethambutol	Myambutol	500,000
Ethamivan		500,000
Ethinamate		500,000

The following compounds tested NEGATIVE on the Methadone 300 ng/mL assay.

Negative Compounds	Trade Name	Concentration Tested (ng/mL)
Ethinyl Estradiol	Estinyl	500,000
Famotidine	Pepcid, Pepcid AC	500,000
Fenfluramine HCl	Pondimin	500,000
Fenopropfen	Nalfon	500,000
Fentanyl citrate	Duragesic	500,000
Flufenamic acid		500,000
Flunitrazepam	Rohypnol	500,000
Fluoxetine	Prozac, Sarafen	500,000
Fluphenazine HCl	Prolixin	500,000
Flurazepam	Dalmane, Durapam	500,000
Flurbiprofen	Ansaid	500,000
Furosemide	Lasix	500,000
Gentamicin	Garamycin	500,000
Gentisic acid		500,000
Glutethimide		500,000
Glycobenclamide		500,000
Glycopyrrolate	Robinul	500,000
Griseofulvin	Fulvicin, Grisactin, Gris-PEG	500,000
Guaiacol glycerol ether (Guaifenesin)	Actifed, Drixoral, Robitussin, Hytuss	500,000
Guanethidine	Ismelin	500,000
Halazepam	Paxipam	500,000
Haloperidol	Haldol, Peridol	500,000
Heroin HCl		500,000
Hexamethylenetetramine		500,000
Hexobarbital		500,000
Hippuric acid		500,000
Hydralazine	Apresoline	500,000
Hydrochlorothiazide	Esidrix, Hydro-DIURIL, Hydro-D, Oretic	500,000
Hydrocodone bitartrate	Lortab	500,000
Hydromorphone HCl	Dilaudid	500,000
Hydroxycoumarin		500,000
Hydroxyphenobarbital		500,000
Hydroxyprogesterone		500,000
Hydroxyzine 2HCl	Anxanil, Apo-Hydroxyzine, Atarax, Hyzine-50	500,000

The following compounds tested NEGATIVE on the Methadone 300 ng/mL assay.

Negative Compounds	Trade Name	Concentration Tested (ng/mL)
Ibuprofen	Advil, Haltran, Mediprene, Motrin, Nuprin, Proen, Rufen	500,000
Idapamine		500,000
Iminostilbene		500,000
Imipramine	Tofranil	500,000
Indapamide	Lozol	500,000
Indole-2-carboxylic acid		500,000
Indole-3-acetic acid	Heteroauxin, Rhizopin	500,000
Indole-3-butyric acid	Hormodin, Seradix	500,000
Indomethacine	Novomethacin, Indomethacin	500,000
Insulin chain $\beta$ , oxidized	Humalog, Humuline, Iletin, Lantus, Novolin, NovoRapid, NovoLog, Oralin, Velosulin	500,000
Ipratropium Br	Attrovent	500,000
Iproniazid		500,000
Isoetharine	Bronkomrter, Bronkosol, Dey-Dose, Dey-Lute, Isoetharine	500,000
Isoniazid	Laniazid, Tubizid, INH	500,000
Isoproterenol HCl	Isuprel	500,000
Isoxsuprine HCl	Vasodilan, Voxsuprine	500,000
Kanamycin	Kantrex	500,000
Ketamine HCl	Ketalar	500,000
Ketoprofen	Orudis, Oruvail	500,000
Labetalol	Normodyne, Trandate	500,000
Lamotrigine	Lamictal	1,000,000
LAMPA		500,000
Levorphanol tartrate	Levo-Dromoran	500,000
Levothyroxine (T <sub>4</sub> )	Synthroid	50,000
Lidocaine	Xylocaine	500,000
Lisinopril	Zestril	500,000
Lithium carbonate		500,000
Lorazepam	Ativan	500,000
Lormetazepam		500,000
Loxapine succinate	Loxitane	500,000
Lysergic acid		500,000
Lysergic Acid Diethylamide (LSD)		500,000
Maprotiline HCl	Ludiomil	500,000

The following compounds tested NEGATIVE on the Methadone 300 ng/mL assay.

Negative Compounds	Trade Name	Concentration Tested (ng/mL)
Mebendazole	Mebendacin, Mebutar, Nemasole, Vermox	500,000
Medazepam	Nobrium	500,000
Medroxyprogesterone acetate	Amen, Curretab, Cycrin, Depo-Provera, Provera	500,000
Mefenamic acid	Ponstel, Ponstan	500,000
Melphalan	Alkeran	500,000
Meperidine	Demerol	500,000
Mephentermine hemisulfate	Wyamine	250,000
Mephénytoin	Mesamtoin	500,000
Mephobarbital	Mebaral	500,000
Mepivacaine HCl	Carbocaine, Isocaine, Polocaine	500,000
Mescaline	Peyote	500,000
Metanephrine HCl		500,000
Metapyriline		500,000
Metaproterenol	Alupent, Metaprel	500,000
Metaraminol bitartrate	Aramine	500,000
Methamphetamine	Desoxyn	500,000
Methaqualone HCl	Normi-Nox, Pallidan, Somnomed, Quaalude	500,000
Metharbital	Gemonil	500,000
Methenamine	Hiprex, Mandameth, Mandelamine, Urex	500,000
Methocarbamol	Robamol, Delaxin	500,000
Methohexital	Brevital Sodium	500,000
Methotrimeprazine	Levoprone	500,000
Methoxamine HCl	Vasoxyl	500,000
Methoxyphenamine HCl		500,000
Methsuximide	Celontin	500,000
Methyldopa	Aldomet	500,000
Methylphenidate	Ritalin	500,000
Methypylon	Noludar	500,000
Metoclopramide	Reglan, Octamide	500,000
Metoprolol tartrate	Lopressor	500,000
Metronidazole	Flagyl, Metizol, Protostat	500,000
Miconazole	Monistat	500,000
Midazolam	Versed	500,000

The following compounds tested NEGATIVE on the Methadone 300 ng/mL assay.

Negative Compounds	Trade Name	Concentration Tested (ng/mL)
Monoethylglycinexylidide HCl (MEGX)		500,000
Morphine		100,000
Morphine-3-Glucuronide		100,000
Morphine-3-β D-Glucuronide		100,000
Nadolol	Corgard	500,000
Nafcillin	Unipen	500,000
Nalidixic acid	Neggram	500,000
Nalorphine HCl		500,000
Naloxone HCl	Narcan	500,000
Naltrexone HCl	Depade, ReVia	500,000
Naphazoline HCl	AK-Con, Albalon, Clear Eyes, Comfort, Degest 2, Estivin II, I-Naphline, Muro's Opcon, Nafazair, Naphcon, Naphcon Forte, Napholine, Ocu-Zoline, Opcon, VasoClear	500,000
Naproxen	Aleve, Anaprox, Naprosyn	500,000
Nelopam HCL		500,000
Neomycin	Mycirfradin Sulfate, Neo-fradin, Neo-Tabs, Myciguent	500,000
Niacinamide (Nicotinamide)		500,000
Nialamide		500,000
Nicotine	Habitrol, Nicoderm, Nicorette, Nicotrol, ProStep	500,000
Nicotinic acid (Niacin)		500,000
Nifedipine	Procardia	500,000
Nitrazepam		500,000
Nitrofurantion	Marcobid	500,000
Nitrophenol		500,000
Norchlordiazepoxide		500,000
Norcodeine HCl		500,000
Nordiazepam (Desmethyldiazepam)		500,000
Nordouripramine		500,000
Nordoxepin HCl		500,000
Norepinephrine HCl	Levophed	500,000
Norethindrone	Aygestom, Micronor, Norlutin	500,000
Norfloxacin	Noroxin	500,000
Norfluoxetine HCl		1,000,000

The following compounds tested NEGATIVE on the Methadone 300 ng/mL assay.

Negative Compounds	Trade Name	Concentration Tested (ng/mL)
Norgestrel	Ovrette	500,000
Normorphine		500,000
Noroxymorphone		500,000
Norpropoxyphene		500,000
Nortriptyline HCl	Aventyl, Pamelor	500,000
Noscapine HCl		500,000
Nylidrin HCl		500,000
Octopamine HCl		500,000
Ofloxacin	Floxin	500,000
Omeprazole	Prilosec	1,000,000
Orotic Acid		500,000
Orphenadrine HCl	Flexoject, Myophen, Norflex	500,000
Oxazepam	Serax	500,000
Oxycodone HCl	OxyContin	500,000
Oxymetazoline HCl	Dristan, Sinarest	500,000
Oxymorphone	Numorphan	500,000
Oxyphenbutazone		500,000
Pargyline	Eutron	500,000
Penicillin G	Pfizerpen	500,000
Penicillin V potassium	PenVee K, Robicillin VK, Veetids, V-Cillin K	500,000
Pentazocine HCl	Talwin	500,000
Pentobarbital	Nembutal	500,000
Pentoxifylline	Trental	500,000
Perphenazine	Trilafon	500,000
Phenacetin		500,000
Phenazopyridine HCl	Azo-Standard, Pyridium, Urodine	500,000
Phencyclidine (PCP)		500,000
Phendimetrazine	Bontril, Plegine	500,000
Phenelzine	Nardil, Parnate, Marplan	500,000
β-Phenethylamine HCl		500,000
Phenformin		500,000
Pheniramine	Trimeton	100,000
Phenobarbital	Luminal	500,000
Phenothiazine		500,000
Phenoxazine		500,000

The following compounds tested NEGATIVE on the Methadone 300 ng/mL assay.

Negative Compounds	Trade Name	Concentration Tested (ng/mL)
Phenoxy methyl penicillic		500,000
Phentermine	Adipex, Fastin, Ionamin, Obenix, Obephen, Obermine, Obermine, Obestin, Phentamine	500,000
Phensuximide		500,000
Phenylalamine		500,000
Phenylbutazone	Azorid, Butazolidin	500,000
Phenylcyclohexylamine		500,000
Phenylcyclopropylamine		500,000
Phenylenediamine		500,000
Phenylpropanolamine	Dura-Vent, Acutrim, Dexatrim	500,000
Phenyltoloxamine citrate salt	Poly-Histine-D	500,000
Phenytoin		500,000
Phthalic acid		500,000
Picrotoxin		500,000
Pipemidic acid		500,000
Piroxicam	Feldene	500,000
Potassium chloride	Kay Ciel, Micro-K, Kaochlor, K-Dur, K-Lor, Kaon-CL, Slow-K, K-Norm, Ten-K	500,000
Potassium iodide	Pima, SSKI	500,000
Prazepam	Centrax	500,000
Prazosin HCl	Minipress	500,000
Prednisolone	Hydeltrasol (Injection)	500,000
Prednisone	Deltasone, Meticorten, Orasone, Stempred	500,000
Pregnenolone		500,000
Prilocaine HCl	Citanest	500,000
Primidone	Mysoline	500,000
Probenecid	Benemid, Benuryl, Probalan	500,000
Procainamide HCl	Procan SR, Procamide SR, Promine, Pronestyl, Pronestyl-SR, Procanbid	500,000
Procaine HCl	Novocain	500,000
Procyclidine	Kemadrin	250,000
Progesterone	Crinone, Ovrette, Micronor, Nor-QD	500,000
Promazine	Primazine, Protine-50, Sparine	100,000
Promethazine HCl	Phenergan	500,000
Propoxyphene	Darvon	500,000

The following compounds tested NEGATIVE on the Methadone 300 ng/mL assay.

Negative Compounds	Trade Name	Concentration Tested (ng/mL)
Propranolol	Betachron, Inderal	500,000
Propylbenzene		500,000
Protriptyline HCl	Vivactil	500,000
Pseudoephedrine HCl	Sudafed	500,000
Pyridinium chlorochromate		500,000
Pyridoxine (Vitamine B <sub>6</sub> )	Nestrex	500,000
Pyrilamine maleate salt		500,000
Quinidine	Cardioquin, Cin-Quin, Duraquin, Novoquinidin, Quinaglute Dura-Tabs, Quinalan, Quinatime, Quinidex Extentabs, Quinora, Quin-Release	500,000
Quinine HCl	Quinamm, Quindan, Quiphile, Q-Vel, Strema	500,000
Quinine sulfate	Formula Q	500,000
Ranitidine	Zantac	500,000
Rescinnamine	Moderil	500,000
Reserpine	Novoresperine, Resertia, Serpalan, Serparil	500,000
Retinoic acid		500,000
Rifampicin	Rifadin	500,000
Ritodrine HCl	Yutopar	500,000
Salicylic acid	Disalcid, Mono-Gesic, Salflex, Salgesic, Salsitab	500,000
Salicyluric acid		500,000
Scopolamine	Transderm-Scop	500,000
Secobarbital	Seconal	500,000
Spirolactone	Novospiroton, Aldactone	1,000,000
Strychnine		500,000
Sulfamethazine Na <sup>+</sup>		500,000
Sulfamethoxazole	Bactrim, Cotrim, Septra, Sulfamethoprim, Sulfatrim, Uroplus	500,000
Sulfathiazole	Sulfa-Gyn, Sulnac, Sultrin, Trysul, V.V.S.	500,000
Sulfisoxazole	Gantrisin	500,000
Sulindac	Clinoril	500,000
Talbutal		500,000
Temazepam	Restoril, Razepam, Temaz	500,000
Terbutaline	Brethaine, Brethine, Bricanyl	500,000

The following compounds tested NEGATIVE on the Methadone 300 ng/mL assay.

Negative Compounds	Trade Name	Concentration Tested (ng/mL)
Terfenadine	Seldane	500,000
Tetracaine HCl	Ak-T-Caine, CepacolViractin, Pontocaine, Opticaine	500,000
Tetracycline	Achromycin V, Parimycin, Sumycin, Tetracyn, Robitet	500,000
Tetrahydrocortisone		500,000
Tetrahydrozoline HCl	Tyzine, Visine	500,000
Thanzylamine		500,000
Theophylline	Aerolate, Asmalix, Elixophyllin, Quinbron, Respbid, Slo-Phyllin, Slo-bid, Theo-24, Theo-Dur, Theolair	500,000
Thiamine HCl	Betalins, Biamine	500,000
Thiothixene	Navane	500,000
Thyroxine		500,000
Timolol	Betimol, Blocadren, Timoptic	500,000
Tocopherol phosphate		500,000
Tolbutamide	Orinase	500,000
Tolmetin	Tolectin	500,000
Toluic acid		500,000
Tramadol	Ultram	500,000
Tranlycypromine HCl	Parnate	500,000
Trazodone HCl	Desyrel	500,000
Triamcinolone	Aristocort, Aristospan, Azmacort, Kenalog, Nasacort	500,000
Triamterene	Dyrenium	500,000
Triazolam	Halcion	500,000
Tribenzylamine		500,000
Triflupromazine HCl	Vesprin	200,000
Trihexylphenidyl HCl	Artane	500,000
Trimethoprim	Prloprim, Trimplex	500,000
Trimipramine	Surmontil	500,000
Tripelennamine HCl	PBZ, PBZ-SR	500,000
Tripolidine HCl	Actagen-C, Allerfrin w/Codenine, Aprodine w/Codenine, Triacin-C, Trifed-C	500,000
Tropic acid		500,000
Tropine	l-Tropine, Ocu-Tropine	500,000
Tryptamine HCl		500,000

The following compounds tested **NEGATIVE** on the Methadone 300 ng/mL assay.

Negative Compounds	Trade Name	Concentration Tested (ng/mL)
Tryptophan		500,000
Tyramine		500,000
Tyrosine		500,000
Urea	Aquacare, Carmel, Grodon's Urea, Ureaphil	500,000
Uric acid		500,000
Venlafaxine	Effexor	1,000,000
Verapamil	Isoptin, Verelan	500,000
Warfarin	Coumadin	500,000
Zolpidem hemitartrate	Ambien	50,000
Zomepirac		500,000
Zopiclone	Imovane	500,000